



2826K

> Perimeter Workstations

Provides work space for up to four-students while storing your laboratory equipment in these perimeter workstations. Choose from 4 different configurations (see below). The four cabinets which make up a single unit measure 24"W x 16"D x 35"H. The stations are designed to utilize space along walls where they can stand alone or be ganged together. Durable solid oak and oak veneer construction finished with a chemical resistant, earth-friendly UV finish. Each unit includes: your choice of 1" solid epoxy resin

or 1" solid phenolic resin top with back splash, solid epoxy drop-in style sink with strainer, 2 combination cold water and gas fixtures, 4 rod sockets, rod storage (rods sold separately), 2 GFI electrical receptacles, keyed-alike locks on all doors and dovetailed drawers and rubber base molding. Assembly and installation of electricity and plumbing is required. Sink measures: 12"W x 8"D x 6"H. Unit measures 90"W x 84"D x 36"H. Work surface on wall portion is 18"D.

ITEM #	W/SINK & FIXTURES	TOP	CABINET STYLE
2814K	Yes	Solid Phenolic	Door
2824K	Yes	Solid Phenolic	Door/Drawer Unit
2834K	Yes	Solid Phenolic	Door/4 Drawer Unit
2844K	Yes	Solid Phenolic	4 Drawer Unit
2816K	Yes	Solid Epoxy	Door
2826K	Yes	Solid Epoxy	Door/Drawer Unit
2836K	Yes	Solid Epoxy	Door/4 Drawer Unit
2846K	Yes	Solid Epoxy	4 Drawer Unit
2814KF	No	Solid Phenolic	Door
2824KF	No	Solid Phenolic	Door/Drawer Unit
2834KF	No	Solid Phenolic	Door/4 Drawer Unit
2844KF	No	Solid Phenolic	4 Drawer Unit
2816KF	No	Solid Epoxy	Door
2826KF	No	Solid Epoxy	Door/Drawer Unit
2836KF	No	Solid Epoxy	Door/4 Drawer Unit
2846KF	No	Solid Epoxy	4 Drawer Unit
100604* / Optional Upright Rod Set w/Full Crossbar			

Note: Units are built to be installed against a wall and do not have finished backs.



2816K



2836K



2846K



2944K | Front View



2944K | Back View



> Forward Vision II

Everything science teachers want. Designed with teacher input and desire for a safe, secure classroom. Teachers maintain control of the class with all students facing the instructor. The large work surface—46"D x 88"W—is made of chemical resistant phenolic resin and holds the solid epoxy sink, and two combination cold water & gas fixtures. Data and electrical outlets are provided on the side of the unit. The cabinets provide locked storage space. Side mounted hooks hold book bags off the floor. The units are finished with a chemical resistant, earth-friendly UV finish. Includes two GFI AC duplex outlets. Also available with a flat top (no sink or fixtures). This model is compact enough to fit multiple units in a small classroom and is easy to assemble. Dimensions 88"W x 46"D x 36"H. Wt. 503 lbs.

STYLE	ITEM #	SINK DIMENSIONS	TOP
Forward Vision II	2944K*	12" x 8" x 6"	Solid Phenolic
Forward Vision II w/Flat Top	2944KF	N/A	Solid Phenolic



Pedestal Base Option

> Octagon Workstation

Crafted of solid hardwood/oak veneers, and finished with a chemical resistant, earth-friendly UV finish. Both units provide ample space for 4 people and are available in 56" or 62", solid 1" epoxy resin or phenolic resin tops. The pedestal with apron comes equipped with 4 drawers, the drawer cabinet unit comes equipped with 16 drawers. Both units come with 4 GFI AC duplex receptacles, and can be ordered with these options: two combination cold water/gas fixtures, 4 rod sockets and, a 16" x 16" x 7" solid epoxy sink with strainer, stopper and trap. The 62"W units with drawers have 2 pulls per drawer. Easily accessible to the utilities. All units are 36"H. (Shipped KD)

Inside Drawer Dimensions:

- 1516K/1514K includes 4 each (20-7/16"W x 10-3/4"D x 2-23/32"H)
- 1546K/1544K includes 16 each (20-7/16"W x 10-3/4"D x 4-15/32"H)
- 1616K/1614K includes 4 each (26-7/16"W x 10-3/4"D x 2-23/32"H)
- 1646K/1644K includes 16 each (26-7/16"W x 10-3/4"D x 4-15/32"H)



Drawer Base Option

WITH SINK	FLAT TOP	STYLE	SIZE	TOP
1516K	1516KF	Pedestal	56"W	Solid Epoxy
1546K	1546KF	Drawer	56"W	Solid Epoxy
1616K	1616KF	Pedestal	62"W	Solid Epoxy
1646K	1646KF	Drawer	62"W	Solid Epoxy
1514K	1514KF	Pedestal	56"W	Solid Phenolic
1544K	1544KF	Drawer	56"W	Solid Phenolic
1614K	1614KF	Pedestal	62"W	Solid Phenolic
1644K	1644KF	Drawer	62"W	Solid Phenolic

> Ez-Lift Workstation (Manual & Push-Button)

Both of the units' height is adjustable from 28" to 40"H. Choose from 1" epoxy resin or phenolic resin top. Students have plenty of space with either the 96"W x 50"D or 88"W x 46"D work surface. Two hooks mounted underneath hold up to four book bags. The grey powder-coated steel and oak veneer base conceals the lifting mechanism. It is finished with a chemical resistant, earth-friendly UV finish. The unit provides six GFI protected electrical outlets. Assembly and power source required. The manual- crank unit has a fold-away crank accessible on the student side of the workstation. The push-button unit has an easy, power option. The motor is accessible by removing one of the oak panels. Plug only supplies power to the hydraulics on the push-button model. AC Duplex needs to be hard-wired in the field.



Push-Button

ITEM #	TOP SIZE	TOP	ADJUSTABLE UNIT
2994K	96"W x 50"D	Solid Phenolic	Manual-Crank
2996K	96"W x 50"D	Solid Epoxy	Manual-Crank
2904K	88"W x 46"D	Solid Phenolic	Manual-Crank
2906K	88"W x 46"D	Solid Epoxy	Manual-Crank
2974K	96"W x 50"D	Solid Phenolic	Push-Button
2976K	96"W x 50"D	Solid Epoxy	Push-Button
2984K	88"W x 46"D	Solid Phenolic	Push-Button
2986K	88"W x 46"D	Solid Epoxy	Push-Button



Fold-Away Crank



> Ez-Lift Workstation

This flexible unit is designed to be used by four students, all facing the instructor. The fold-away crank allows the units height to be adjusted from 28"H to 40"H, for ease of use with a wheelchair. The front of the unit has locking storage measuring 36"W x 8"D x 24"H. Inside the cabinet are two lipped shelves and a quad electrical outlet. This space can be used for computer storage, recharging or for general storage. The front panel includes two more duplex outlets and two data ports. The unit is finished with a chemical resistant, earth-friendly UV finish and boosts 500 lb. weight capacity. The top is 1" solid phenolic resin. Overall size: 96"W x 54"D x 28"-40"H. Assembly and power source required.

ITEM #	SIZE	TOP	ADJUSTABLE UNIT
2894K	96"W x 54"D x 28" - 40"H	Solid Phenolic	Manual-Crank



C2716K



C2516K



(C2726K)

Two-student island workstations are fitted together to form a continuous row of workstations.

> Two-to-Eight Student Island Workstations

Single faced tables that provide a fully equipped science center from two to eight students. Each unit is constructed of solid oak, and hardwood veneers. The single faced units contain a solid epoxy resin sink (12" x 8" x 6"H) along with cold water faucets for every two students. Every unit contains enough gas cocks for each student along with GFI AC duplex electrical receptacles, rod sockets and book compartments. All units contain locks for the drawers and doors. Choose from models with a full cupboard only or a half cupboard and dovetailed drawers. Can be ordered without sink, fixtures and rod sockets. They are finished with a chemical resistant, earth-friendly UV finish. Single face two-student unit measures 66"L x 30"W x 36"H. Add 66" in length for each additional two students. Your choice of 1" solid epoxy resin or 1" phenolic resin tops.

W/SINK FIXTURES & ROD SOCKETS	W/O SINK FIXTURES & ROD SOCKETS	STUDENTS #	TOP	CABINET STYLE
C2514K	C2514KF	2	Solid Phenolic	Drawers
C2524K	C2524KF	4	Solid Phenolic	Drawers
C2534K	C2534KF	6	Solid Phenolic	Drawers
C2544K	C2544KF	8	Solid Phenolic	Drawers
C2516K	C2516KF	2	Solid Epoxy	Drawers
C2526K	C2526KF	4	Solid Epoxy	Drawers
C2536K	C2536KF	6	Solid Epoxy	Drawers
C2546K	C2546KF	8	Solid Epoxy	Drawers
C2714K	C2714KF	2	Solid Phenolic	Cupboard
C2724K	C2724KF	4	Solid Phenolic	Cupboard
C2734K	C2734KF	6	Solid Phenolic	Cupboard
C2744K	C2744KF	8	Solid Phenolic	Cupboard
C2716K	C2716KF	2	Solid Epoxy	Cupboard
C2726K	C2726KF	4	Solid Epoxy	Cupboard
C2736K	C2736KF	6	Solid Epoxy	Cupboard
C2746K	C2746KF	8	Solid Epoxy	Cupboard



C2416K



C2616K



(C2426K)

Two-student island workstations are fitted together to form a continuous row of workstations.

> Four-to-Sixteen Student Island Workstation

Double faced tables that provide a fully equipped science center from four to sixteen students. Each unit is constructed of solid oak and hardwood veneers and finished with a chemical resistant, earth-friendly UV finish. The double faced units contain a solid epoxy resin sink (12" x 8" x 6"H) along with cold water faucets for every two students. Every unit contains enough gas cocks for each student along with GFI AC duplex electrical receptacles, rod sockets and book compartments. All units contain locks for the drawers and doors. Choose from models with a full cupboard only or a half cupboard and dovetailed drawers. Can be ordered without sink, fixtures and rod sockets. Double face four student unit measures 66" L x 48" W x 36" H. Add 66" in length for each four students.

W/SINK FIXTURES & ROD SOCKETS	W/O SINK FIXTURES & ROD SOCKETS	STUDENTS #	TOP	CABINET STYLE
C2414K	C2414KF	4	Solid Phenolic	Drawers
C2424K	C2424KF	8	Solid Phenolic	Drawers
C2434K	C2434KF	12	Solid Phenolic	Drawers
C2444K	C2444KF	16	Solid Phenolic	Drawers
C2416K	C2416KF	4	Solid Epoxy	Drawers
C2426K	C2426KF	8	Solid Epoxy	Drawers
C2436K	C2436KF	12	Solid Epoxy	Drawers
C2446K	C2446KF	16	Solid Epoxy	Drawers
C2614K	C2614KF	4	Solid Phenolic	Cupboard
C2624K	C2624KF	8	Solid Phenolic	Cupboard
C2634K	C2634KF	12	Solid Phenolic	Cupboard
C2644K	C2644KF	16	Solid Phenolic	Cupboard
C2616K	C2616KF	4	Solid Epoxy	Cupboard
C2626K	C2626KF	8	Solid Epoxy	Cupboard
C2636K	C2636KF	12	Solid Epoxy	Cupboard
C2646K	C2646KF	16	Solid Epoxy	Cupboard