

MOBILETECH

LAB TABLES

The **MOBILE TECH LAB TABLE** is fabricated using high performance materials for sustainability. Flexibly designed and carefully constructed to provide a durable, streamlined solution for the modern lab, the LS-Mobile Tech Lab Table eliminates the need for other bulky lab furniture. It's the perfect workstation for any sustainable lab environment where flexibility is mandatory for routine work efficiency.

FEATURES

- + **Powder Coated Steel Frame** [inset 3-6" or flush frame design]
- + **5" Heavy Duty Locking Casters**
- + **HPL // Chemical Resistant Top** [high pressure laminate]
- + **Work Surface Sizes** [2', 2.5', 4', 5' or 6']
- + **Welded Steel Frame** [1200lb. weight capacity]
- + **Worktop Cutout for Cable Management** [concealment of wires + hoses]



MOBILETECH

A STREAMLINED
SOLUTION FOR
THE MODERN LAB.



Designed for workspace adaptability, Mobile Tech Tables allow easier access to all services and provide more usable space in any lab environment. Select from a variety of frame colors, durable work surfaces and table sizes. Accessories such as under-bench drawer storage and adjustable shelving are also available.

MOBILE TECH TABLE
2' Worktop



MOBILE TECH TABLE
2.5' Worktop



MOBILE TECH TABLE
5' Worktop



MOBILE TECH TABLE
6' Worktop



MOBILE TECH TABLE
4' Worktop

