

# FIRM FLEX

## LAB CHAIRS

The **FIRM FLEX LAB CHAIR SERIES** is puncture-proof, tear resistant and designed to last. Durable constructed of the same material as a modern car dashboard, the LS-Firm Flex Chair also provides firm lumbar support, comfort where it counts and adjustability options to ensure ergonomic correctness for each user. Flex lab chair features include independent seat back, height + tilt pressure levers, as well as a fixed or mobile base option.



Adjustable Height from 17-25"

Nylon, Chrome or Aluminum Base Style

Fixed or Mobile - 5 Leg Base

Puncture-Proof Polyurethane Seat



# FLEXCHAIRS

## FEATURES

- + **Height Adjustable** [seat back + tilt pressure on 3 lever models]
- + **Fixed or Mobile Base Style** [nylon, chrome or aluminum]
- + **5-Leg Standard Base** [glides or 2.5" casters]
- + **Durable Seat Material** [tear resistant Polyurethane]
- + **Heavy Duty Welding** [holds up to 450lb.]

## SEAT SURFACES



**POLYURETHANE COLORS**  
BLACK

## BASE STYLES

### MOBILE



**NYLON 25"W**



**CHROME 26"W**



**ALUMINUM 27"W**

### FIXED



**NYLON 25"W**



**CHROME 26"W**



**ALUMINUM 27"W**

## DIMENSIONS

ITEM #	SEAT	HEIGHT RANGE	BACK	BASE	BASE STYLE
<b>LS-LNT-UA</b>	18"D X 18.5"W	18 to 26"H	12"H X 16"W	25"W	Nylon
<b>LS-LCT-UA</b>	18"D X 18.5"W	21" to 29"H	12"H X 16"W	26"W	Chrome
<b>LS-LAT-UA</b>	18"D X 18.5"W	17" to 25"H	12"H X 16"W	27"W	Aluminum
<b>LS-LNT-UB</b>	18"D X 18.5"W	22" to 32"H	12"H X 16"W	25"W	Nylon
<b>LS-LCT-UB</b>	18"D X 18.5"W	25" to 35"H	12"H X 16"W	26"W	Chrome
<b>LS-LAT-UB</b>	18"D X 18.5"W	22" to 32"H	12"H X 16"W	27"W	Aluminum