

LAB TABLES

METAL FLEX SERIES

The **FLEX SERIES LAB TABLE** has the versatility to accommodate all disciplines including General Science, Technology, Robotics, Engineering and Math. Flexibly designed to meet the unique requirements of the modern lab, each table is durably constructed with industrial grade, heavily welded steel.



Height Adjustable 30-36"

Heavy-Duty Steel Frame

Fixed or Mobile Base Style

Easy Access to Services

FLEXTABLES

FEATURES

- + **Height Adjustable** [30-36"// sit to stand table position]
- + **Fixed Style or Mobile Base** [4" heavy duty locking casters]
- + **Powder Coated Steel Frame** [available in a variety of colors]
- + **Durable Work Surface Options** [phenolic resin or plastic laminate]
- + **Service Supply Modules** [easy access to H/C water, gas, electricity & data]
- + **Under-Bench Storage** [optional drawer, cupboard or combination storage]
- + **Adjustable Reagent Shelving** [with optional LED task lighting]



TABLE FRAMES



STEEL COLORS [60% Gloss]

WHITE | BLACK | STORM GRAY | BEIGE | BLUE | DARK BLUE



WORK SURFACES



PHENOLIC RESIN COLORS

TITAN WHITE | BLACK



PLASTIC LAMINATE COLORS

WHITE | BLACK | CLOUDY GRAY | BEIGE | STORM GRAY

DIMENSIONS

ITEM #	TABLETOP	HEIGHT ADJUST	FRAME	WORK SURFACE OPTIONS
LS-KY2424	24"D X 24"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2430	24"D X 30"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2436	24"D X 36"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2448	24"D X 48"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2460	24"D X 60"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2472	24"D X 72"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2496	24"D X 96"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate

▶ **ALSO AVAILABLE** 30" | 36" | 48"D